Tool ID: 1808 Tool Location: 224

### **Equipment Information Sheet**

# Laurell WS-400A spin coater (2nd Floor)

Calls to staff phones will be automatically forwarded to Manager: Xinwei Wu 607-254-4934 Backup: Michael Skvarla 607-254-4674

their cell phones during accessible hours. At other times

leave a message or send them an email.

#### **SAFETY**

• splash hazard: wear proper protective gear for working in solvent hoods.

#### **USAGE RESTRICTIONS**

- ALL materials must be cleared FIRST with the tool manager(s).
- Never clean the spinner. Tell the tool owner each day that you use the spinner for cleaning.

#### SCHEDULING/SIGN-UP RESTRICTIONS

Minimum Tool Time: 0 minutes

• None (no need to log into CORAL).

## MATERIALS COMPATIBILITY CATEGORY

Tool Category 5: Class A and B Metals and Compounds	
Allowed	Not Allowed
Tool category 1/1E, 2, 3, and 4 materials	
Silicon Based Substrates and Films	
III/V compound Semiconductors	
Glass Substrates	
PECVD and ALD Films	
Cured organics and baked Photoresist	
CNF Class A, B, and Refractory metals	
Exposed Gold, Silver, Copper	
Alkali and Alkaline Compounds	
Organic/Biology Molecules prepared- w/salt buffers	
High Vapor Pressure Materials (Mg, Ca, Zn)*	* Some tool restrictions on high vapor pressure materials may apply
Soft organic materials	

#### High Vapor Pressure Metals and Compounds are materials that have a vapor pressure above 1e-6 Torr at 400 C.

#### **Additional Material Restrictions and Exceptions**

• Wear eye protection, face shield and thin nitrile gloves when using the spinner in the fume hood.

Last Updated: 01/13/2022