



- The System uses pressurized EBR 10A
- This solvent is pressurized at about 10 PSI
- This setup utilizes a Nordson EBR stainless steel needle valve with it's corresponding controller.
- The solvent dispensing works along with a Brewer Science flange spinner model 100.



- The cantilever style arm can be adjusted in both lateral and horizontal dimension.
- Only the lateral dimension has a ruler.
- Users should always check with Lithography staff before making any adjustments.



- The cantilever arm can swing away.
- This is only used for spinner maintenance.



- The tip controller is always on and setup.
- After a power shutdown specially during the monthly power test the controller will display this message.



- If you press SEL button once the tool will automatically load tip 1.
- This is the setup that talk to the Brewer spinner micro-controller.



- 18 second is the time the unit will open the valve to release EBR 10A.
- For now this will be the permanent setup for this tool.



- This is the normal display of the spinner's micro-controller.
- The LCD will display
 - Run Process
 - Edit Process
 - Diagnostics
- Pressing Edit Process will the display the recipe that I have setup for EBR




- The EBR recipe consist of two process steps.
- Step 1
 - 300 RPM [Speed], 400 RPM/Sec [Acceleration], 20 Sec [Step Time], #2 [Dispense]
 - Note that Dispense time will be 18 seconds
- Step 2
 - 2000 RPM [Speed], 1400 RPM [Acceleration], 15 Sec [Step Time]



- To load the EBR recipe go back to the original display or just select "Load" from the previous screen.



- To load the EBR recipe go back to the original display or just select “Load” from the previous screen.



- Select the EBR recipe and click “OK”




- The controller will give you a warning asking you if you really want to load the recipe.




- Press “YES”



- This screen will be displayed.
- It will display “Spin Process EBR”




- Before you load your substrate make sure to clean the excess solvent on the tip using a clean-room swab.
- Since the system is constantly under low pressure it will built up a small EBR 10A liquid ball.



- The lateral dimension has been setup by the CNF staff for 100mm substrates.
- Very carefully load your substrate and center it using the black centering tool.
- Remove centering tool without touching your substrate.




- Press “Center” and ensure that your substrate is centered before running the recipe.
- The black button on the side can also be used to select center.
- If it looks good press the black button again and the controller will run the two steps in the recipe.
- Remember that the Nordson controller will active during the first step of the recipe for 18 second.
- The second step will dry the substrate.



- The edge bead removal will start right after centering.
- Keep in mind that if you press the selector button too fast the spinner controller will skip the spin part of the recipe, in other words it will skip step 1 and step 2.
- If step 1 and step 2 are skipped the EBR 10A will still dispense and ruin your coated substrate.
- Please be careful.



- The EBR recipe can be run multiple times in order to clear thick or old resist.



- Be very careful when removing your substrate.
- Make sure you don't touch the tip.