

Tool ID: 428
Tool Location: 121

Equipment Information Sheet
UV Curing Station

Manager:	Roberto Panepucci	607-254-4855	Calls to staff phones will be automatically forwarded to their cell phones during accessible hours. At other times leave a message or send them an email.
Backup:	Giovanni Sartorello	607-254-4853	

SAFETY

- Do not Lift UV Curing Box while ON!
- Turn it OFF to take your sample out!
- Use handle to lift after finished. It can be HOT!

USAGE RESTRICTIONS

SCHEDULING/SIGN-UP RESTRICTIONS

MATERIALS COMPATIBILITY CATEGORY

Minimum Tool Time: 0 minutes

Additional Material Restrictions and Exceptions

Last Updated: 02/13/2024